



Single Use Anesthesia Breathing Masks

Single Use Anesthesia Breathing Masks from GE Healthcare are tested and verified to function with GE Healthcare anesthesia machines. Anesthesia breathing masks are designed to channel gas mixtures to the patient for non-invasive ventilation via the mouth and nose. The products perform optimally with GE Healthcare anesthesia machine breathing systems when applicable.

FEATURES

- Clear face piece allows visual observation of patient status
- Compatible with all standard breathing systems with 22mm and 15mm connectors where applicable
- Available in sizes from neonatal through large adults
- Single use masks offer pre-filled soft inflatable cushion that can be inflated and deflated to enhance the seal to suit the patient's anatomical contours.
- Single use masks include color coded hooks to help indicate mask size

**Our equipment.
Our accessories.
Simply compatible.**

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

SINGLE PATIENT USE ANESTHESIA MASK OFFERINGS

Pre-filled inflatable single patient use anesthesia masks are designed to channel anesthesia gases or fresh gas mixtures to the patient for inhalation via the mouth and nose. Anesthesia masks are used in conjunction with anesthesia and ventilator breathing systems.

PRODUCT SPECIFICATIONS

Material composition

Component	Material
Cushion	Polyvinyl Chloride (PVC)
Mask Shell	Polyvinyl Chloride (PVC)

Latex content

Pre-Filled Inflatable Disposable Anesthesia Face Masks do not contain natural rubber latex.

DEHP content

Pre-Filled Inflatable Disposable Anesthesia Face Masks do not contain phthalate DEHP.

Product use

Pre-Filled Inflatable Disposable Anesthesia Face Masks are single use devices only.

Sterility

Pre-Filled Inflatable Disposable Anesthesia Face Masks are supplied non-sterile.

Storage

Store in a cool, dry place out of direct sunlight

PART LIST

Part Number	Description	Quantity
2106570-012	Mask, disposable, size 0, 15mm connection	50
2106570-013	Mask, disposable, size 1, 15mm connection	50
2106570-014	Mask, disposable, size 2, 22mm connection	50
2106570-015	Mask, disposable, size 3, 22mm connection,	50
2106570-016	Mask, disposable, size 4, 22mm connection	50
2106570-017	Mask, disposable, size 5, 22mm connection	50
2106570-018	Mask, disposable, size 6, 22mm connection,	50

Shelf life

Shelf life of 5 years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

Disposal considerations

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Packaging materials

Primary	Polybag – Polyethylene (PE)
Secondary	Carton / Box – Cardboard

Product may not be available in all countries and regions.

Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com

Distributed by:
Datex-Ohmeda, Inc.



Manufactured by:
Flexicare Medical Limited
Cynon Valley Business Park,
Mountain Ash, CF45 4ER, UK