




AVURI[®]

One Centralized Solution for Device Management

Remotely manage a fleet of medical devices across
multiple locations with AVURI by GE Healthcare.



EMPOWER YOUR HEALTHCARE TECHNOLOGY MANAGERS

AVURI is a cloud-based solution that lets you manage devices remotely from one central location. AVURI doesn't just allow you to monitor devices across your fleet but it also enables your Healthcare Technology Management teams the ability to view, back up and deploy device configurations, all without leaving their desks.

Maximize Productivity



Remote management can help reduce on-site visits and can help ensure minimum equipment downtime.

Drive Consistency



Standardized clinical configurations and hospital protocols across your medical devices.

Increase Availability



Increase device uptime and improve patient flow.

SAVE TIME, OPTIMIZE RESOURCES

Significant time is currently spent performing tasks device by device that could be improved by using AVURI Device Management.



**3 – 4
hrs/device/yr**

Average time spent performing tasks.



>70%

Must be performed locally on device.

GE Internal Research*

Secure

Every aspect of AVURI is built with security in mind. Your device data on the AVURI platform is protected by GE Healthcare's proprietary medical grade data security standards and data that is transferred is never patient data.

User Management

AVURI allows central users to assign and control access to others and all actions are logged providing improved visibility.

Customized Views

AVURI is designed to allow you to prioritize actions, view progress and isolate errors. With a single click you can customize your view to match your needs.

Simple Device Setup

Simply activate the feature on your supported GE Healthcare medical device and you are ready to go. With AVURI's self-service device claiming feature, devices are displayed real-time within the AVURI platform.

Local Acceptance

AVURI is designed to not interfere with daily usage. Updates from the cloud to a device first need to be accepted by the local user in order to ensure minimal disruption.

Intuitive User Interface

AVURI displays a unified list of all assets on a single screen. The intuitive UI uses tiles to represent devices, drag and drop functionality and an intuitive sorting and grouping functionality.

* Fleet Manager Opportunity Assessment Quantitative Customer Market Research, 2019.

Product may not be available in all countries and regions. Full product technical specification is available upon request. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com/promotional-locations. Data subject to change.

© 2020 General Electric Company.

GE, the GE Monogram and AVURI are trademarks of General Electric Company.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

September 2020
JB00027XU

