


MotionFree – 1st ever digital respiratory motion management solution

Seamlessly integrate respiratory motion management into your existing workflow



Always on and seamless

No external device
No manual setup



No uncomfortable gating device

No need to coach patients or explain an external system



Routine within existing workflow

Auto-correction only in anatomies that need it



MotionFree delivers up to a 30% improvement in quantitative accuracy (SUV_{mean}) and up to a 67% improvement in volumetric accuracy of lesion size*

*Compared to non-processed (STATIC, no motion correction) data, as demonstrated in phantom testing using a typical and fast respiratory model, 18 mm Ge-68 spheres and OSEM reconstruction.