



SENOGRAPHE 700/800T WITH INTEGRATED SHIELD



MINIMUM:

10'-0" x 12'-0" [3.0m x 3.65m]

7'-8" [2.3m]



9-18

① 1 SENDGRAPHE 700/800T MAMMOGRAPHY GANTRY WITH BASIC SHIELD AND CONSOLE 665 lbs 4952 btu

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

- ⑥0 MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS
 24 IN. W x 77 IN. H [610mm x 1956mm], CONTINGENT ON
 A 60 IN. [1524mm] CORRIDOR WIDTH
- ⑥1 COUNTER TOP WITH BASE AND WALL CABINETS
- ⑨0 MAIN DISCONNECT CONTROL, WALL MOUNTED, GEMS CAT. NO.
 E4502EE, 65 lbs., SEE DETAIL R4502EE.
- ⑨1 EDM TOWER FOR SEND 700/800T

50-75° F (10-24° C)
 15° F (10° C) per hour
 30-75 percent non condensing
 10 percent per hour
 9500 ft. (29,000 M) above sea level