


Vascular Care with a Personal View

# OEC Outpatient Vascular Premium Suite



## Why own an Outpatient Vascular suite?

- Control your schedule and your time
- Know and care for your patients personally
- Grow your practice and your brand



Only GE can offer this premium turnkey solution to help you transition to a more personal level of vascular care...for you and your patients.

- A suite of premium products that operate like a fixed interventional suite
- A world-class dedicated support team
- Financing options including low monthly payments



# Design your Premium Suite


## Image Visualization and Integration



- Displays multi-modality images on up to 6 Hi-Def monitors
- Positionable monitors on articulating arms for comfortable viewing
- Fully installed in 2 days

NuBOOM Premium Visualization and digital communication system

## C-arm Choices



**Performance: 9800 Plus Vascular C-arm**

- 12" large field of view
- 15 kW high-powered generator for penetration
- Rotating anode X-ray tube for case throughput
- DSA & Roadmap for vascular imaging



**Premium: 9900 Elite Vascular & Motorized C-arms**

- 12" large field of view
- 15 kW high-powered generator for penetration
- Rotating anode X-ray tube for case throughput
- DSA & Roadmap for vascular imaging



- DRM/MTS designed for vascular imaging
- Motorized C-arm moves like a fixed room system
- Table side control (RUI) to control C-arm functions
- Active tube cooling for additional throughput

## Table Choices



**IDI - Aspect ISR**

- 4-Way Floating Tabletop
- Isocentric Lateral Tilt ("Iso-Roll" Patent Pending)



**Stille - imageQ2**

- 90" table top
- 4 way floating table top
- 660 lbs. patient weight cap.
- Exchangeable table top

## Include an Injector



**MedRad Injector Arterion**

- Latest technologies
- Light, maneuverable
- Easy to use

Finance your entire Suite with GEHC Financial Services to enjoy one manageable monthly payment.



imagination at work

©2012 General Electric Company - All rights reserved.

GE Healthcare,  
Surgery - Americas:  
Phone 801-328-9300  
Fax 801-328-4300

LR-050571-01