

Ultrasound

U.S. Accessories & Supplies



How to Order



Shop online:
gehealthcare.com/accessories

Available 24/7, GE Healthcare Service Shop makes ordering supplies online fast and easy.

- Search real-time inventory of compatible, OEM products
- Take advantage of Auto Reorder for parts you buy regularly
- Batch order multiple items at once with Quick Order
- Get instant order and tracking confirmations
- Pay by PO or credit card or integrate your purchasing system

Visit GE Healthcare Service Shop to register.

1. Visit GE Healthcare Service Shop, navigating to the registration link in the upper left-hand corner of the page.
2. Give us your account details. Pro tip: Provide your customer account number for expedited registration.
3. A welcome message will come to you through the email you provide. Start shopping!

Register now at ServiceShop.GEHealthcare.com



Other ways to order:

Phone: (866) 281-7545 Options 1, 2, 1, 2

Fax: (877) 279-6990

Email: gehaccessorysales@ge.com

EDI: Free – contact us for details

Hours: 7AM – 6PM (CST)



Table of Contents

Transducer Disinfection.....	U/S 4	▶
Transducer Storage.....	U/S 8	▶
trophon Products and Accessories.....	U/S 11	▶
Endocavity/General Probe Storage & Covers.....	U/S 14	▶
ABUS System Consumables.....	U/S 16	▶
SensorySuite for ABUS.....	U/S 17	▶
General Purpose Transducer Covers.....	U/S 20	▶
Endocavity Transducer Covers.....	U/S 22	▶
Ultrasound Lotion, Disinfecting Wipes & Gels.....	U/S 24	▶
Endocavity Needle Guides.....	U/S 27	▶
Furniture.....	U/S 29	▶
Printer Paper & Accessories.....	U/S 30	▶
Power Quality.....	U/S 31	▶

TD 200 portfolio of products to clean and disinfect TEE probes

- **Fastest automated cycle time in the US:** 10 minute automated disinfectant time for TEE probes; with total cycle time improved from 17 minutes to 10 minutes
- **Fastest High Level Disinfection of TEE Probes cleared by the FDA:** 3 minutes; high level disinfection time improved from 5 minutes to 3 minutes
- Single Use TD-12 AquaCide disinfectant – no longer an Aldehyde based solution – therefore a neutralization station is not needed



TD 200 Automated TEE Probe Cleaner Disinfectant



E8364GA



E8364GB



E8364GC



E8364GF



E8364GE



E8364GD

Cat#	Description	Warranty
E8364GA	TD 200 Automated TEE Probe Cleaner Disinfectant	C
E8364GB	TD-12 AquaCide	W
E8364GC	QwikCheck® Chemical Indicators	W
E8364GF	TD 200 HLD Tracker™ Documentation Log	W
E8364GE	AquaCide Spill Kit	W
E8364GD	Replacement vapor management and drain filter – TD 200	W

TD 100 portfolio of products to disinfect and protect TEE probes

TD 100 TEE Probe High Level Disinfecter

The TD 100® Automated TEE Probe Disinfecter is designed to provide high-level disinfection of transesophageal (TEE) echocardiogram ultrasound probes.

- Simple to use: quick 5 min disinfection cycle, in addition to 5 automated rinse cycles for total of 17 min
- Closed system with advanced vapor management: captures and neutralizes disinfectant fumes
- The automatic waste disposal feature helps to eliminate the potential for over-soaking the probe
- Repeatable, verified disinfection cycle: printed record for proof of successful cycle
- Tested and approved by the major TEE manufacturers



TD 100 TEE Probe Disinfecter



E8364BA



E8364BG



E8364DH



E8364DJ



E8364BE

Cat#	Description	Warranty
E8364BA	TD 100 Disinfecter (preinstallation kit, starter kit, vapor management filter & one case of disinfectant)	C
E8364BG	TD 100 Thermal Paper Roll (each)	W
E8364BF	Replacement Vapor Management Filter for TD 100 and Drain Filter	W
E8364DH	Nephros Water Filter Starter Kit	W
E8364DJ	Nephros Annual Maintenance Kit (includes 2 filters)	W
E8364BE	TD-8 Single-use, 0.59% ortho-pht halaldehydebased HLD, 32/box	W

TEEClean portfolio of product to clean and disinfect TEE probes

TEEClean Automated TEE Probe Cleaner Disinfector

TEEClean® is an automated solution that both cleans and high-level disinfects TEE probes. Automated reporting – including electronic record data logging and verification reporting – facilitates regulatory compliance/audits.

- Eliminates the manual cleaning process and need for a sink
- Reduces the potential for probe damage due to oversoaking
- Minimizes probe handling with integrated electrical leak testing
- Extra step of MRC testing NOT required – and single use disinfectant means no exposure to an open disinfectant that must be changed out regularly



E8365FM

TEEClean Automated TEE Probe both cleans and high-level disinfects. It eliminates the manual cleaning process and need for a sink.

TEEClean Automated TEE Probe Cleaner Disinfector



E8365FM



E8364FN



E8364FP



E8364FQ



E8364BE

Cat#	Description	Warranty
E8365FM	TEEClean Automated TEE Probe Cleaner Disinfector	C
E8364FN	TEEZyme for TEEClean	W
E8364FP	Nephros 5 mm filter (includes 4 filters)	W
E8364FQ	Vapor Management Filtration System	W
E8364BE	TD-8 Single-use, 0.59% ortho-pht halaldehydebased HLD, 32/box	W



TEEClean, TD 100 & TD 200 Quick Guide

Cat#	Description	Use with TD 200	Use with TD 100	Use with TEEClean
E8364GA	TD 200 Automated TEE Probe Disinfector	X		
E8364BA	TD 100 Disinfector		X	
E8364FM	TEEClean Cleaner & Disinfector			X
E8364GB	TD-12 AquaCide	X		
E8364BE	TD-8 Single-use Cartridge		X	X
E8364GC	QwikCheck® Chemical Indicators	X		
E8364FN	TEEZyme for TEEClean			X
E8364FF	TEEZyme®MC Enzymatic Cleaner	X	X	
E8364DJ	Nephros Annual Maintenance Kit	X	X	
E8364GD	Replacement Vapor Management and Drain Filter – TD 200	X		
E8364BF	TD 100 Vapor Management and Drain Filter		X	
E8364FP	Replacement Water Filter 5mm			X
E8364FQ	Main and Drain Vapor Management filter - TEEClean			X
E8364BM	TEEZyme Contoured Enzymatic Sponge	X	X	X
E8364BP	TEEZyme Sponge Storage Unit	X	X	X
E8364DF	QwikDry® Probe Drying Cloth	X	X	X
E8364DG	QwikDry Storage Unit	X	X	X
E8364DE	Electrical Leak Tester	X	X	X
E8364BQ	Leakage Tester Lead Wire	X	X	X
E8364DS	Printer Kit for Leakage Tester	X	X	X
E8364DT	Wall Mount for Leakage Tester & Printer	X	X	X
E8364EC	Printer Paper for Leakage Tester	X	X	X
E8364BC	Neutralization Station		X	X
E8364BJ	Neutralization Station Replacement Carbon Filter		X	X
E8364BH	Glycinex™		X	X
E8364GE	AquaCide Spill Kit	X		
E8364BN	Master Spill Kit		X	X
E8364GF	TD 200 HLD Tracker™ Documentation Log	X		
E8364EF	TD 100 HLD Tracker™ Documentation Log		X	
E8364BG	Thermal Paper Roll	X	X	X
E8364BB	CleanShield® TEE Probe Storage Cabinet (3 probe)	X	X	X
E8364BD	CleanShield® TEE Probe Storage Cabinet (6 probe)	X	X	X
E8364CC	HEPA Replacement Filter for Storage Cabinet	X	X	X
E8364CD	PreFilter Replacement Pack for Storage Cabinet	X	X	X
E8364DV	TPorter® Probe Transportation & Procedure Case	X	X	X
E8364EH	TPorter Cart	X	X	X
E8364FE	TPorter Handle Carrier	X	X	X
E8364DZ	PullUp™ TEE Probe Bio-Barrier Sleeve (red)	X	X	X
E8364DY	PullUp TEE Probe Cover (clear)	X	X	X
E8364DX	TPorter Tie Bands	X	X	X

Probe Care and Storage - TEEClean, TD 100 or TD 200



E8364BB



E8364BD



E8364CC



E8364CD



E8364BM



E8364BP



E8364DF



E8364DG

Cat#	Description	Warranty
E8364BB	CleanShield® TEE Probe Storage Cabinet (3 probe capacity) 88 in (H) x 25.25 in (W) x 10 in (D), 115 lb	C
E8364BD	CleanShield TEE Probe Storage Cabinet (6 probe capacity) 88 in (H) x 50.5 in (W) x 10 in (D), 230 lb	C
E8364CC	HEPA Replacement Filter for CleanShield Ultrasound Storage Cabinet (1 filter/box)	W
E8364CD	PreFilter Replacement Pack for CleanShield Ultrasound Storage Cabinet (12 filter/box)	W
E8364BM	TEEZyme Contoured Enzymatic Sponge for TEE probes, Box of 100 Sponges	W
E8364BP	TEEZyme Sponge Storage Unit	W
E8364DF	QwikDry Probe Drying Cloth, Case of 100 Cloths	W
E8364DG	QwikDry Storage Unit	W

Compliance Products - TEEClean, TD 100 or TD 200



E8364DE



E8364DT



E8364DS



E8364BN



E8364BQ



Probe adapter for leakage tester.
Contact us for details

Cat#	Description	Warranty
E8364DE	Electrical Leak Tester Base Model with added real time features (requires lead and adapter)	C
E8364BQ	Leakage Tester Lead Wire	C
E8364DT	Wall Mount for Leakage Tester & Printer	C
E8364DS	Printer Kit for Leakage Tester	C
E8364EC	Printer Paper for Leakage Tester	W
E8364BN	Master Spill Kit	W

Transport - TEEClean, TD 100 or TD 200



E8364DV



E8364EH



E8364FE



E8364DY



E8364DZ



E8364DX

Cat#	Description	Warranty
E8364DV	TPorter® Probe Transportation and Procedure Case	C
E8364EH	TPorter Cart Holds 4 Tporter cases for transport	C
E8364FE	TPorter Carrier Carrier hooks Tporter case to console handle	C
E8364DY	PullUp TEE Probe Cover 100 sleeves/box	W
E8364DZ	PullUp™ TEE Probe Bio-Barrier Sleeve 100 sleeves/box	W
E8364DX	TPorter Tie Bands Pack of 200 (100-red, 100-green)	W

Neutralization - TD 100, TEEClean



E8364BC



E8364BH

Cat#	Description	Warranty
E8364BC	Neutralization Station	C
E8364BH	Glycinex, 2 oz packets (24 per case)	W
E8364BJ	Neutralization Station Replacement Carbon Filter	W

trophon²[®] - Automated High Level Disinfecter

trophon² is an environmentally friendly, efficient system for the complete High-Level Disinfection (HLD) of ultrasound transducers. It features AcuTrace™ (RFID) technology, which enhances digital record keeping to deliver a cost-effective, easy and reliable workflow for your HLD process.

trophon Sonex-HL™

Catalog #E8350MC

The cartridge contains a proprietary solution of 35% hydrogen peroxide. 80 mL each. Product cannot ship air.

trophon Chemical Indicators

Catalog #E8350MB

Chemical Indicators confirm that there is a Minimum Effective Concentration of disinfectant.



E8350MA



E8350MC



E8350MB

Cat#	Description	Packaging	Warranty
E8350MA	trophon2 Device	X	P
E8350MB	trophon Chemical Indicator (1 in/26 mm)	300/package	P
E8350MC	trophon Sonex-HL	6 bottles/box (80 ml bottle)	P

Traceability Accessories

trophon AcuTrace Operator Card

Catalog #E8350MJ

Incorporates AcuTrace™ technology to digitally capture operator interaction details, stored directly on your trophon2.



E8350MJ

trophon AcuTrace Medical Instrument Tag

Catalog #E8350MK

Incorporates AcuTrace™ technology to digitally capture ultrasound probe details, stored directly on your trophon2.



E8350MK

trophon AcuTrace PLUS License

Catalog #E8350ML

AcuTrace™ PLUS allows you to create a customized solution for streamlined workflows.



E8350ML

Cleaning and Drying Wipes

trophon Drying Wipes

Catalog #E8350MR

Are used to dry probes before and after HLD. The Drying Wipes are highly absorbent to leave surfaces completely dry.



E8350MR

trophon Wipes Wall Mount

Catalog #E8350MS

For wall mounting trophon Cleaning or Drying Wipes.

trophon Cleaning Wipes

Catalog #E8350MQ

The Wipes provide fast and efficient low level disinfection (LLD) of ultrasound equipment.



E8350MQ

Probe Storage Solutions

Clean Ultrasound Probe Covers

Catalog #E8350TA

Protect endocavitary and surface ultrasound probes from the risk of recontamination.



E8350TA

Cat#	Description	Packaging	Warranty
E8350MJ	trophon AcuTrace Operator Card	10/box	P
E8350MK	trophon AcuTrace Medical Instrument Tag	10/box	P
E8350ML	trophon AcuTrace PLUS License	1	P
E8350MQ	trophon Companion Cleaning Wipes	6/box (100 wipes/can)	P
E8350MR	trophon Companion Drying Wipes	6/box (100 wipes/can)	P
E8350MS	trophon Companion Wipes Wall Mount	1	P
E8350TA	Clean Ultrasound Probe Cover	100 covers/box	P

trophon Mounting Options

trophon2 Wall Mount2

Catalog #E8350MD
 12.0 in (L) x 11.5 in (W) x 10.8 in (H)
 Weight – 8.4 Lb/3.8 kg

trophon Mobile Cart

Catalog #E8350ME
 36.4 in (H) x 20.7 in (W) x 22.4 in (D)
 Weight – 53 Lb/24 kg



E8350MD



E8350ME

trophon Printer and Accessories

trophon Printer

Catalog #E8350PA
 trophon Printer Product Dimension:
 8.23 in (L) x 4.29 in (W) x 6.73 in (H)

Records each HLD cycle and generates an adhesive label for documentation and traceability.

trophon Printer Paper/Labels (1 roll – 600 prints)

Catalog #E8350PB
 Labels are direct thermal adhesive with a seven year print stability.



E8350PA



E8350PB

Printer Mounting Options

trophon Printer Wall Mount

Catalog #E8350PF
 Custom designed for secure, horizontal mounting of the trophon Printer including cable organization.

Printer Mount for trophon Cart 2

Catalog #E8350PG
 Designed for securely attaching the trophon Printer to the trophon Cart.
 (Note: only for customers who have both the trophon Cart and Printer).



E8350PF



E8350MG

Cat#	Description	Packaging	Warranty
E8350MD	trophon Wall Mount 2	1	P
E8350ME	trophon Cart	1	P
E8350PA	trophon Printer	1	P
E8350PB	trophon Printer Paper	1 roll – approx 600 prints	P
E8350PF	trophon Printer Wall Mount	1	P
E8350MG	Printer Mount for trophon Cart 2	1	P

Endocavity/General Purpose Probe Storage

CleanShield® Probe Storage Cabinet

Designed to effectively and securely store disinfected endocavity probes and semi-critical devices.

- Specifically designed for trans vaginal, trans rectal, and general purpose ultrasound probes
- Stores up to 6 endocavity probes
- Probes stored in vertical position to minimize stress on connection cable, strain relief and electrical pack
- Made of thermally-fused polypropylene – easy to clean and will not rust
- Unique two-stage filtration removes potential airborne contaminants from entering the cabinet
- Dimensions – 88 in (H) x 26 in (W) x 10 in (D)
- Weight – 115 lbs.



E8364CA



E8364CA

CleanShield® Ultrasound Storage Cabinet - HEPA Filter

- Recommended filter replacement interval is once every 2 years, or when indicator light comes on
- 1 filter/box



E8364CC



E8364CD

CleanShield® Ultrasound Storage Cabinet, PreFilter pack

- Recommended filter replacement interval is once every 3 months
- 12 filters/box

Cat#	Description	Packaging	Warranty
E8364CA	CleanShield® Ultrasound Storage Cabinet, 110V	1	C
E8364CC	HEPA Filter for CleanShield® Ultrasound Storage Cabinet	1/box	W
E8364CD	PreFilter Pack for CleanShield® Ultrasound Storage Cabinet	12/box	W

Ultrasound Probe Rack

Our Ultrasound Probe Rack has been designed to hold a variety of probes.

- Cable support on back plate can hold small, narrow probes ie: MSK
- Includes 4 pre-drilled holes for wall mounting
- Weight – 13.5 lbs
- Measurements – 20 in (L) x 15 in (H) x 5.5 in (D)



E8363JF

Cat#	Description	Warranty
E8363JF	Ultrasound Probe Rack	H

Compatibility: Long intracavity, large 4-D OB and standard size probes

Disposable Probe Storage Covers

Protek Non-sterile Probe Covers

- Fits endocavity and general purpose probes
- Helps protect entire probe, including handle
- Designated space for applying disinfection label
- Latex-free



E8388AA



E8388AB

Cat#	Description	Packaging	Warranty
E8388AA	Protek Non-Sterile Disposable Storage Cover – 4" x 16" clear storage cover	100/box 4 boxes/case	H
E8388AB	Protek Non-Sterile PullUp™ Disposable Storage Cover – 4" x 16" PullUp Storage cover with Secure Tape	200/box	H

Clean Non-sterile Ultrasound Probe Covers

Decrease Risk of Recontamination after High Level Disinfecting Probes Clean Ultrasound Probe Covers protect intracavity and surface ultrasound probes from the risk of recontamination through handling or the environment after high level disinfection (HLD).

- Enclosed wire wrap fastens cover around the probe and handle
- Internal cover size 4.3 in (W) x 13.6 in (L)
- Manufactured in a certified clean room (ISO 14644-1, Class 8)
- Latex-free
- Single-use
- Non-sterile
- 100/box



E8350TA

Cat#	Description	Packaging	Warranty
E8350TA	Clean Ultrasound Probe Covers	100/box	P

ABUS System Consumables

ABUS Membranes

The ABUS disposable scan membrane provides comfort to the patient during the 3D image acquisition. The flexible, polyester, mesh fabric stabilizes the breast tissue and enables better image quality during the scan. The membrane is single-use-only which eliminates the risk of cross-contamination and the need for sterilization.

ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS and somo•v™ ABUS Platinum

Catalog #E8340AB

Curved membrane dimensions:

229 mm (L) x 213 mm (W) x 19 mm (D)

Polysonic® Ultrasound Lotion and Protex™ Ultrasound Disinfectant Wipes

See details on page U/S 18



E8340AB

GE Healthcare offers an array of consumable supplies for the ABUS (Automated Breast Ultrasound) portfolio of systems including scan membranes, Polysonic® ultrasound lotions, and Protex™ disinfectants.



E8340AB

ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo•v™ ABUS (Automated Breast Ultrasound) Platinum systems (non-Sterile)



E8365EF

E8365EC

Compatible for use with either ABUS membrane



E8365EG

E8365EB

Compatible for use with ABUS system

Cat#	Description	Quantity per Box	Warranty
E8340AB	ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo•v™ ABUS (Automated Breast Ultrasound) Platinum systems (non-Sterile)	50/box	W
E8365EC	Polysonic Ultrasound Lotion (Non-Sterile) – 8.5 oz dispenser bottles	12 bottles/box	H
E8365EF	Polysonic Ultrasound Lotion (Non-Sterile) – 1 US gallon with refillable dispenser	4 gallons/case	H
E8365EB	Protex Ultra Disinfectant Wipes – 60 ct Soft-pack w/stay-put adhesive backing – 6.5 in x 6.0 in	12 pack/box	H
E8365EG	Protex Ultra Disinfectant Wipes – 75 ct canister pack – 7 in x 9.5 in	8 pack/case	H

SensorySuite for ABUS

An all encompassing interactive experience centered around your patient.



Sight & Sound

A flat-screen monitor will display soothing video and relaxing ambient sounds of your selected environment.

Sight

Decorative wall panels in a choice of three ambiances and two sizes help set the mood for the selected ambiance.

Sight

A flat-screen monitor on the wall will display a slide show of photos of your selected ambiance.

Scent

A scent diffuser gently infuses the air with SensorySuite's light calming fragrance.

By focusing on your patients, we've re-invented the experience... and addressed many of the reasons why women avoid breast imaging.

Designed to stimulate a woman's senses, SensorySuite can distract her from the perceived discomfort, pain, and anxiety of a breast exam.

SensorySuite lets your patient choose the environmental ambiance she prefers for her breast imaging:

- Seaside
- Garden
- Waterfall
- Mountain

Designed to help you stand apart, SensorySuite helps enable you to:

- Improve the breast imaging experience of every woman
- Increase patient compliance with screening breast imaging
- Promote the high level of patient care you provide
- Differentiate your healthcare facility from others in your area

Advertise and promote your SensorySuite solution to women and referring physicians in your area with our comprehensive marketing tool kit of promotional materials.

Take breast imaging to a whole new place!



Seaside



Garden



Waterfall

SensorySuite may be purchased as packages or as individual items.



Seaside



Waterfall



Garden

Platinum Package includes:

- Interactive Experience including two large monitors (diag. 48 in/122 cm)
- SensorySuite Scent diffuser
- SensorySuite Exam Room fragrance capsules (Set of 30 scent diffuser refill capsules)
- One Wall Panel – Large size (106 cm x 49 cm/41.1 in x 19.3 in)
- Three Wall Panels – Large size 3 x (106 cm x 49 cm/41.1 in x 19.3 in)

Cat#	Description	Warranty
E6301RA	SensorySuite – Platinum Package – Seaside	E
E6301PA	SensorySuite – Platinum Package – Waterfall	E
E6301SA	SensorySuite – Platinum Package – Garden	E

Silver Package includes:

- Interactive Experience including one large monitor (diag. 48 in/122 cm)
- One Wall Panel – Large size (106 cm x 49 cm/41.1 in x 19.3 in)

Cat#	Description	Warranty
E6301SC	SensorySuite – Silver Package – Seaside	E
E6301RC	SensorySuite – Silver Package – Waterfall	E
E6301PC	SensorySuite – Silver Package – Garden	E

Gold Package includes:

- Interactive Experience including two small monitors (diag 32 in/81 cm)
- SensorySuite Scent diffuser
- SensorySuite Exam Room fragrance capsules (Set of 30 scent diffuser refill capsules)
- Three Wall Panels – Large size 3 x (106 cm x 49 cm/41.1 in x 19.3 in)

Cat#	Description	Warranty
E6301PB	SensorySuite – Gold Package – Seaside	E
E6301RB	SensorySuite – Gold Package – Waterfall	E
E6301SB	SensorySuite – Gold Package – Garden	E

Bronze Package includes:

- Interactive Experience including one small monitor (diag 32 in/81 cm)
- One Wall Panel – Small size (70 x 32 cm/27.5 x 12.6 in)

Cat#	Description	Warranty
E6301RD	SensorySuite – Bronze Package – Seaside	E
E6301PD	SensorySuite – Bronze Package – Waterfall	E
E6301SD	SensorySuite – Bronze Package – Garden	E

SensorySuite may be configured with individual items.

SensorySuite - Sight - Patient interactive device

Cat#	Description	Warranty
E63301CE	Tablet PC for waiting room	Y
E63301CF	Secured and flexible wall mount system for tablet PC in the waiting room	W



E63301CE

SensorySuite - Sight - Decorative Wall Panels

Cat#	Description	Warranty
E63301AA	Waterfall - One Wall Panel - Large size (106 x 49 cm/41.1 x 19.3 in)	W
E63301AB	Waterfall - Three Wall Panels - Large size 3 x (106 x 49 cm/41.1 x 19.3 in)	W
E63301AC	Waterfall - One Wall Panel - Small size (70 x 32 cm/27.5 x 12.6 in)	W
E63301AD	Waterfall - Three Wall Panels - Small size 3 x (70 x 32 cm/27.5 x 12.6 in)	W
E63301AE	Seaside - One Wall Panel - Large size (106 x 49 cm/41.1 x 19.3 in)	W
E63301AF	Seaside - Three Wall Panels - Large size 3 x (106 x 49 cm/41.1 x 19.3 in)	W
E63301AG	Seaside - One Wall Panel - Small size (70 x 32 cm/27.5 x 12.6 in)	W
E63301AH	Seaside - Three Wall Panels - Small 3 x (70 x 32 cm/27.5 x 12.6 in)	W
E63301AJ	Garden - One Wall Panel - Large size (106 x 49 cm/41.1 x 19.3 in)	W
E63301AK	Garden - Three Wall Panels - Large size 3 x (106 x 49 cm/41.1 x 19.3 in)	W
E63301AL	Garden - One Wall Panel - Small size (70 x 32 cm/27.5 x 12.6 in)	W
E63301AM	Garden - Three Wall Panels - Small size 3 x (70 x 32 cm/27.5 x 12.6 in)	W



SensorySuite - Scent: Diffuser/Capsules

Cat#	Description	Warranty
E63301BA	SensorySuite Scent diffuser	Y
E63301BB	SensorySuite Exam Room fragrance capsules (Set of 30 scent diffuser refill capsules)	W



E63301BA



E63301BB

SensorySuite - Sight & Sound - Interactive Experience

Cat#	Description	Warranty
E63301CJ	Interactive Experience including one small monitor (diag 32 in/81 cm)	E
E63301CK	Interactive Experience including two small monitors (diag 32 in/81 cm)	E
E63301CL	Interactive Experience including one large monitor (diag 48 in/122 cm)	E
E63301CM	Interactive Experience including two large monitors (diag 48 in/122 cm)	E

General Purpose Transducer Covers

PullUp™ Sterile Latex-Free General Purpose Transducer Covers

- Cover helps to provide protection against cross-contamination and maintains a sterile field
- The innovative design of the PullUp allows for the cover to be easily opened, easy application of gel, and the easy application of the cover
- Available in lengths of 24 and 48 in
- Packaged 24 Covers per box, sold per box



E8324PC/E8348PC

Cat#	Description	Packaging	Warranty
E8324PC	24 in PullUp Sterile Latex-Free Probe Cover 5.5 × 24 in (14 × 61 cm)	24/box	K
E8348PC	48 in PullUp Sterile Latex-Free Probe Cover 5.5 × 48 in (14 × 122 cm)	24/box	K

Intraoperative Transducer Covers

Protek Sterile Latex-Free Polyethylene Intraoperative Transducer Covers

- Intraoperative transducer cover for a sterile environment
- Extends up the length of the transducer cord, allowing up to 8 ft (2.4 m) of sterile protection
- Catalog #E8333AR include sterile gel packets and bands
- Packaged 12 Covers per box, sold per box



E8333AG



E8333AR

Cat#	Description	Packaging	Warranty
E8333AG	Sterile 6 in (tapers to) 3 in × 96 in (9.14 × 244 cm) Polyethylene Transducer Cover (No Gel)	12/box	K
E8333AR	Sterile 6 in (tapers to) 3 in × 96 in (9.14 × 244 cm) Polyethylene Transducer Cover (Gel)	12/box	K

Protek Sterile Latex-Free Transducer Covers

- Provides an extended field of sterile protection against cross-contamination
- Form-fitting tip provides a snug fit for small footprint transducers
- Packaged 10 Covers per box, sold per box



E8333AS

Cat#	Description	Packaging	Warranty
E8333AS	96 in PullUp Sterile Latex-Free Probe Cover 5.875 x 96 in (12.7 x 243.8 cm) with Form Fitted Tip and 2 Securing Clips	10/box	K

Compatibility: GE i12SL and i12L-RS transducers

Endocavity Transducer Covers – Non-Sterile

Endocavity Probe Covers – Non-Sterile

- Used for OB/GYN scanning procedures
- Helps to protect against cross-contamination
- Available for Auto-Reorder

Non-Sterile

	Cat#	Latex Free	Description	Bands Included	Gel Included	Packaging
	E8333AP	X	12 in (30.48 cm)	X		50/bag
	E8333KB	X	12 in (30.48 cm) with PullUP	n/a		200/bag
	E8333AL	X	2.5 x 24 in (6.4 x 61 cm)	X		24/box
	E8333DL		0.8 x 7.9 in (2 x 20 cm)	X	X	24/box
	E8333BM		1 x 7.9 in (2.6 x 20 cm)	X	X	24/box
	E8333AE		1.4 x 7.9 in (3.5 x 20 cm)	X		24/box
	E8333BR		1.4 x 7.9 in (3.5 x 20 cm)			100/bag
	E8333BS		1.4 x 7.9 in (3.5 x 20 cm)			500/bag
	E8333DM		0.8 x 2.8 in (2 x 7 cm)	200/bag		100/bag
	E8386CG	X	2.37 in tapering to 1.75 in (W) x 9.5 in (L) (60 mm tapering to 44 mm (W) x 241 mm (L))		X	100/bag

Endocavity Probe Cover Bands

	Cat#	Latex Free	Description	Packaging
	E8333AJ	X	.75 in (1.9 CM)	500/bag

NOTE: All Endocavity Probe Covers have Warranty Code K.

Endocavity Transducer Covers – Sterile

Endocavity Probe Covers – Sterile

- Used for OB/GYN scanning procedures
- Helps to protect against cross-contamination
- Available for Auto-Reorder



	Cat#	Sterile	Latex Free	Description	Bands Included	Gel Included	Packaging
	E8333KA	X	X	12 in (30.48 cm) with PullUP	n/a		200/bag
	E8333AK	X	X	2.5 X 24 in (6.4 X 61 cm)	X	X	24/box
	E8333EC	X		0.8 x 7.9 in (2 x 20 cm)	X		24/bag
	E8333BN	X		1.4 x 7.9 in (3.5 x 20 cm)	X		24/bag

Ultrasound Lotion

Parker Labs Polysonic® Ultrasound Lotion

POLYSONIC Ultrasound Lotion is an acoustic coupling lotion that is hypo-allergenic and has a light pleasing scent.



E8365EF

E8365EC

Cat#	Description	Packaging	Warranty
E8365EC	PolySonic Ultrasound Lotion (Non-Sterile) – 8.5 oz dispenser bottles	12/box	H
E8365EF	PolySonic Ultrasound Lotion (Non-Sterile) – 1 US gallon	4 gal/case	H

Compatibility: Used with GE ABUS system

Ultrasound Disinfectant Wipes

Parker Labs Protex Ultrasound Disinfectant Wipes

Protex™ one-step cleaner disinfectant kills more than 30 pathogens including HIV, H1N1 and MRSA, while remaining safe for even the most sensitive surfaces. The pre-moistened wipes are easy and convenient to use.



E8365EG

E8365EB

Cat#	Description	Packaging	Warranty
E8365EB	Protex Ultra Disinfectant Wipes – 60 ct Soft-pack w/stay-put adhesive backing – 6.5 in x 6.0 in	12 pack/box	H
E8365EG	Protex Ultra Disinfectant Wipes – 75 ct canister pack – 7 in x 9.5 in	8 pack/case	H

Compatibility: Visit gehealthcare.com/ultrasound for detailed information.

Ultrasound Gels

Parker Labs Aquasonic® 100 Gel

- Original gel for diagnostic and therapeutic applications
- Tested and approved for use with GE transducers
- Blue colored aqueous gel formula is non-staining, non-sensitizing, and non-irritating
- Hypoallergenic, bacteriostatic
- No salt or formaldehyde; not a spermicide
- 5 liter jug with dispenser, 250 ml bottles, or 20 g sterile overwrapped foil pouches



E8365BA

E8365AF

E8365BJ

Cat#	Description	Packaging	Warranty
E8365BA	Aquasonic 100 Gel (250 ml bottles)	12/case	H
E8365AF	Aquasonic 100 Gel (5 L jug – includes dispenser)	1/case	H
E8365BJ	Aquasonic 100 Gel (20 g pouches – single-use)	48/box	H

Parker Labs Aquasonic Clear Gel

- Same formula as original Aquasonic 100 Gel, but without the perfume or color
- Clear, non-scented aqueous, water-soluble gel formula is non-staining, non-sensitizing and non-irritating
- Hypoallergenic, bacteriostatic
- Offered in 250 ml bottles or 5 liter jug with dispenser



E8365CB

E8365CA

Cat#	Description	Packaging	Warranty
E8365CA	Aquasonic Clear Gel (250 ml bottles)	12/case	H
E8365CB	Aquasonic Clear Gel (5 L jug – includes dispenser)	1/case	H

Parker Labs SCAN® Gel

- Dual viscosity formula allows SCAN gel to be applied as a gel and scan as a fluid
- Less viscous than Aquasonic 100 for ultrasound procedures where a less viscous gel is desired
- Green colored gel is completely water-soluble, will not stain clothing, damage equipment, or irritate skin
- Hypoallergenic, bacteriostatic
- Offered in 250 ml bottles, 1 gallon jug with dispenser



E8365BC

E8365BD/BK

Cat#	Description	Packaging	Warranty
E8365BC	SCAN Ultrasound Gel (8 oz. bottles)	12/case	H
E8365BD	SCAN Ultrasound Gel (1 gal. jug)	1/case	H

Standoff Pads

Parker Labs Aquaflex® Standoff Pads

- Aqueous, flexible, disposable ultrasound standoff pad for use in difficult to visualize and near field areas
- Requires no gel, except in prolonged procedures
- Cushions sensitive areas for comfortable scanning
- Hypoallergenic, bacteriostatic
- Pads are disposable, which helps to prevent cross-contamination
- Each pad measures 2 × 9 cm
- Packaged 6 pads per box, sold per box



E8317C

Cat#	Description	Packaging	Warranty
E8317C	Aquaflex Standoff Pads – 2 × 9 cm	6/box	H

Endocavity Needle Guides

Protek Disposable Endocavity Needle Guide (Single Piece)

- 24 disposable procedure kits for ultrasound guided biopsy
- Disposable design of guide helps to reduce the risk associated with cross-contamination
- Directs instruments according to on-screen biopsy guidelines
- Needle Guides include (E8333JA):
 - 24 disposable biopsy needle guides
- Each starter kit includes (E8333JB):
 - 24 disposable biopsy needle guides
 - 24 latex transducer covers 1.4 × 7.9 in (3.5 × 20 cm) with gel and bands



E8333JA



E8333JB

Cat#	Description	Packaging	Warranty
E8333JA	E8C/E721/E8C-RS/IC5-9H/MTZ Needle Guides	24/box	K
E8333JB	E8C/E721/E8C-RS/IC5-9H/MTZ Biopsy Kit	24 procedures/box	K

Compatibility: E8C, E721, E8C-RS, IC5-9H and MTZ Transducers

Protek Disposable Endocavity Needle Guide

- 25 disposable procedure kits for ultrasound guided biopsy
- Disposable design of guide reduces risks associated with cross-contamination




E8333DK

Cat#	Description	Packaging	Warranty
E8333DK	Sterile TV/TR 618E Endocavity Probe Long Needle Guide	25 procedures/box	K

Compatibility: TB/TR 618E transducers

Protek Endocavity Needle Guides - Disposable

- Procedure kits for ultrasound guided biopsy
- Directs instruments according to on-screen biopsy guidelines
- Each Disposable Biopsy Needle Guide Starter Kit (E8387M) includes 24 sterile kits:
 - Each sterile kit includes: 1 disposable needle guide, 1 latex probe cover, 1 non-latex probe cover and elastic bands
- Disposable design of guide reduces risks associated with cross-contamination
- Cost-effective disposable kits helps maximize patient throughput



E8387M

Protek Endocavity Needle Guide - Reusable

- Stainless steel reusable biopsy guide allows for repeat usage after sterilization requirements have been met per hospital or clinic protocol
- Directs instruments according to on-screen biopsy guidelines



E8387MA

Cat#	Description	Packaging	Warranty
E8387M	BE9C Disposable Biopsy Needle Guide Starter Kit	24 procedures/box	H
E8387MA	BE9C Reusable Biopsy Needle Guide	1	H

Compatibility: BE9C transducers

Ultrasound Furniture

Ultrasound Carrying Case for Venue* 40

- Dent-resistant, shatter-resistant, lightweight, air and watertight case
- Durable molded-in hasps that allow the addition of standard size padlocks
- Hinges with integrated feet that hold fast under impact, and create a flat surface for upright stability



E8365CS

Cat#	Description	Warranty
E8365CS	Venue 40 Carrying Case	H

GE Anti-Fatigue Mat

- Helps to reduce body stress and fatigue resulting from standing for long periods
- One of the most comfortable, buoyant, anti-fatigue floor mats on the market
- Durable construction; sealed to prevent moisture absorption and facilitate cleanup



E7058A



E7058AB



E7058AC

Cat#	Description	Warranty
E7058A	GE Anti-Fatigue Floor Mat (Blue 3' x 5' x 5/8") – 20 lb	H
E7058AB	GE Anti-Fatigue Floor Mat (Gray 3' x 5' x 5/8") – 20 lb	H
E7058AC	GE Anti-Fatigue Floor Mat (Blue 2' x 3' x 5/8") – 9 lb	H

Physician Mobile Armchair

- Padded arms
- Charcoal gray
- Adjusts from 17.7 to 21.7 in (45 to 55 cm)
- Weight 31 lb (14 kg)
- Chair contains some ferro-magnetic material
- Not safe for use in an MRI Suite's Zone IV, use in control room only



E80331BE

Cat#	Description	Warranty
E80331BE	Physician Mobile Armchair	H

Sony Printer Media Supplies

Sony Printer Media and Supplies

- GE offers a variety of high-quality Sony printer papers and supplies to meet all of your printing needs. Please refer to the attached matrix to find the Sony media recommended for your printer



E8310JT



E8310JS



E8310KP



E8310KS



E8310LA



E8310LB



E8310LC



E8310GA

Cat#	Cross-Ref	Description	Compatible with Sony Printers	Prints Per Pack or Roll	Sold Per	Warranty
E8310JT	UPP-210HD	Sony High-Density Black & White Media (130 prints per roll; 1 case contains 5 rolls)	UP-991AD, UP-971AD, UP-910, UP-930, UP-960, UP-980	130	1 case = 5 rolls	H
E8310JS	UPP-210SE	Sony Standard Black & White Media (130 prints per roll; 1 case contains 5 rolls)	UP-991AD, UP-971AD, UP-910, UP-930, UP-960, UP-980	130	1 case = 5 rolls	H
E8310KP	UPC-21L	Sony A6 Color Print Pack (200 prints of paper and ribbon – 50 sheets x 4 packs)	UP-D25MD, UP-25MD, UP-21MD, UP-D23MD, UP-20	200	1 carton = 6 packs	H
E8310KS	UPC-55	Sony A5 Color Print Pack (2 individual packs of 100 sheets each and 2 ribbons)	UP-D55, UP-D55MD	200	1 carton = 5 packs	H
E8310LA	UPP-110HA	Sony High-Density Black & White Media (240 prints per roll; 1 case contains 10 rolls)	UP-880, UP-890MD, UP-D890	240	1 case = 10 rolls	H
E8310LB	UPP-110HD	Sony High-Density Black & White Media (270 prints per roll; 1 case contains 10 rolls)	UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897MD, UP-D898MD, UP-X898MD	270	1 case = 10 rolls	H
E8310LC	UPP-110HG	Sony High-Glossy Black & White Media (240 prints per roll; 1 case contains 10 rolls)	UP-895MD, UP-D895, UP-897MD, UP-D897MD, UP-D898MD, UP-X898MD	240	1 case = 10 rolls	H
E8310GA	UPP-110S	Sony Standard Black & White Media (270 prints per roll; 1 case contains 10 rolls)	UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897MD, UP-D898MD, UP-X898MD	270	1 case = 10 rolls	H

Power Quality

Uninterruptible Power Supplies (UPS)

- Helps provide reliable, clean, consistent power for Diagnostic Imaging and PACS systems
- Helps to eliminate power anomalies (noise, transients, over and undervoltage), which could damage the system's sensitive computer components
- Enables the DI system to complete scanning after the loss of supply power and allows for the saving of valuable data and orderly system shutdown



E8349PA

Cat#	Description	Warranty
Ultrasound UPS		
E4502M	420 VA Powerware UPS	ZZ
E4502KV	2 kVA UPS	D
E8249PA	1440 VA UPS (Medical Grade)	ZZ

U/S	UPS Cat#	UPS Description
LOGIQ 200	TBD	TBD
LOGIQ 400	E4502KV	2 kVA UPS
LOGIQ 500	E4502KV	2 kVA UPS
LOGIQ L7	E4502KV	2 kVA UPS
LOGIQ 700	E4502KV	2 kVA UPS
LOGIQ P5	E4502M	0.4 kVA UPS (Battery - E4502NA) (10 ft cord - E4502NB)
LOGIQ P6	E4502KV	2 kVA UPS
LOGIQ E9	E8349PA	1440 VA UPS (Medical Grade)
LOGIQ P3	TBD	TBD
LOGIQ S7	TBD	TBD
LOGIQ S8	E8349PA	1440 VA UPS (Medical Grade)
Vivid 3	E4502KV	2 kVA UPS
Vivid 3 Expert	E4502KV	2 kVA UPS
Vivid 3 Pro	E4502KV	2 kVA UPS
Vivid FiVe	E4502KV	2 kVA UPS
Vivid 7	E4502KV	2 kVA UPS
Vivid 7 PRO	E4502KV	2 kVA UPS
Vivid P3	TBD	TBD
Vivid FiVe EchoPAC Review Stations, Options and Accessories	TBD	TBD
Voluson E8/E6	TBD	TBD

Save time. Order accessories online with Service Shop.

Auto-Reorder

Auto-Reorder

Set quantity, schedule frequency

Every 3 Months ▾

Customize the quantity and order frequency that works for you.

 Add to Cart

Add to cart

It's that easy. Your first order has been placed. The product will be auto-delivered at your preferred delivery frequency.



Manage reorders

In advance of each delivery, we will send you a reminder email for your upcoming delivery. Skip your deliveries or cancel your subscriptions at any time.



Quick Order

Know what you need to order? Batch order multiple items at once with Quick Order on Service Shop. Enter a list of item numbers or paste from an excel spreadsheet, select quantity and add to cart. Spend less time ordering!



Payment Options

Pay by PO or credit card or integrate your purchasing system.

Register now at [ServiceShop.GEHealthcare.com](https://www.gehealthcare.com/accessories)

1. Visit GE Healthcare Service Shop, navigating to the registration link in the upper left-hand corner of the page.
2. Give us your account details. Pro tip: Provide your customer account number for expedited registration.
3. A welcome message will come to you through the email you provide. Start shopping!



Specifications

- All specifications are subject to change without notice
- The product and brand names herein are the trademarks or registered trademarks of their respective owners

Pricing

- Please call to obtain current pricing. All pricing and product availability is subject to change without notice
- Price invoiced will be the price in effect at the time of order entry
- Quotations available upon request

Blanket Purchase Orders

You can save time and money by ordering your consumable products with our Blanket Order program. Blanket orders must be scheduled for a year or more and a signed purchase order agreement must be on file. Contact us for additional details.

Restocking Fee

- Returns are subject to, but not limited to, a minimum 15% restocking fee
- Certain items, such as sterile and environmentally controlled products within this catalog, are non-returnable and nonrefundable unless defective and have been noted as such

Product Returns

- Immediately upon receipt, please inspect all deliveries for order accuracy and condition
- For products shipped in error, please obtain return authorization from the Accessories and Supplies Call Center within thirty (30) days of delivery
- Power Quality items are nonreturnable and nonrefundable

Contact Us for General Information and Warranty Terms

In U.S.

GE Healthcare
Accessories and Supplies Customer Support Center
P.O. Box 414, W-544
Milwaukee, WI 53202

Phone: (866) 281-7545 – opt 1, opt 2, opt 1, opt 2

Fax: (877) 279-6990

Online: www.gehealthcare.com/accessories

Email: gehccaccessorysales@ge.com

All other regions

Contact your local GE Representative.



www.gehealthcare.com

Harnessing data and analytics across hardware, software and biotech, GE Healthcare is the \$19 billion healthcare business of GE (NYSE:GE). As a leading provider of medical imaging equipment, with a track record of more than 100 years in the industry and more than 50,000 employees across 100 countries, we transform healthcare by delivering better outcomes for providers and patients.

Imagination at work

GE Healthcare
3000 North Grandview
Waukesha, WI 53188
USA
www.gehealthcare.com

©2021 General Electric Company.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. The results expressed in this document may not be applicable to a particular site or installation and individual results may vary. This document and its contents are provided to you for informational purposes only and do not constitute a representation, warranty or performance guarantee from GE Healthcare.

GE, the GE Monogram, and imagination at work are trademarks of General Electric Company.

GE Healthcare, a division of General Electric Company.

All other product names and logos are trademarks or registered trademarks of their respective companies.

JB00640US

Rev. 01/21