


GE Healthcare

Vascular Care with a Personal View

# OEC Outpatient Vascular Modular Suite



## Why own an Outpatient Vascular suite?

- Control your schedule and your time
- Know and care for your patients personally
- Grow your practice and your brand

Only GE can offer this quality modular turnkey solution to help you transition to a more personal level of vascular care...for you and your patients.


- A suite of premium products-for imaging and viewing
- A world-class dedicated support team
- Financing options including low monthly payments



Visit us on the web at  
[www.gehealthcare.com/outpatientvascular](http://www.gehealthcare.com/outpatientvascular)

# Design your Modular Suite


## Image Visualization and Integration



**IDI MDS**  
**Mobile Display System**

- Mobile Display System (MDS) – portable imaging management system that displays various medical images and pictures
- 26" LCD HD monitors
- Simple video routing and control

## C-arm Choices



**Performance: 9800 Plus**  
**Vascular C-arm**

- 12" large field of view
- 15 kW high-powered generator for penetration
- Rotating anode X-ray tube for case throughput
- DSA & Roadmap for vascular imaging



**Premium: 9900 Elite Vascular & Motorized C-arms**

- 12" large field of view
- 15 kW high-powered generator for penetration
- Rotating anode X-ray tube for case throughput
- DSA & Roadmap for vascular imaging



- DRM/MTS designed for vascular imaging
- Motorized C-arm moves like a fixed room system
- Table side control (RUI) to control C-arm functions
- Active tube cooling for additional throughput

## Table Choices



**IDI - Aspect ISR**

- 4-Way Floating Tabletop
- Isocentric Lateral Tilt ("Iso-Roll" Patent Pending)



**Stille - imageQ2**

- 90" table top
- 4 way floating table top
- 660 lbs. patient weight cap.
- Exchangeable table top

## Include an Injector



**MedRad Injector Arterion**

- Latest technologies
- Light, maneuverable
- Easy to use

**Finance your entire Suite with GEHC Financial Services to enjoy one manageable monthly payment.**



imagination at work

©2012 General Electric Company – All rights reserved.

GE Healthcare,  
Surgery — Americas:  
Phone 801-328-9300  
Fax 801-328-4300

LR-050572-01