



GE Temperature Cables



Easy connection. Accurate measurements.

GE temperature cables support disposable and reusable temperature probes used with GE monitoring systems.

Product compatibility

GE temperature cables are tested and verified with the following GE equipment:

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE Patient Data Module, E-PSMP and E-PT
- B105/B125 patient monitors
- B20/B40 and B20i/ B40i patient monitors

Accessories simply compatible with your equipment.

You can trust GE Healthcare accessories to keep your clinical devices working efficiently.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step



PRODUCT SPECIFICATIONS

Materials

Not made with natural rubber latex, DEHP or BPA
ROHS (Restriction of hazardous substances) compliant

ORDERING INFORMATION

Part Number	Description
2104178-002	Temperature cable, 400 series probes, dual, 0.5 m/1.7 ft. Connects directly to reusable temperature probes.
2104179-001	Temperature cable, 400 series for disposable probes, single, 3.6 m/12 ft. Connects directly to reusable temperature probes.
2104178-001	Temperature cable, 400 series probes, single, 0.5 m/1.7 ft. Connects directly to temperature probes.

Product may not be available in all countries and regions.
Contact a GE Healthcare Representative for more information.
Please visit www.gehealthcare.com/promotional-locations

Imagination at work

© 2018 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram, are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE, the GE monogram, Imagination at work and CARESCAPE are trademarks of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2146744 6/18