



GE Healthcare

A lifetime of reassurance

Corometrics™ Nautilus Nine-Crystal transducers are available directly from GE Healthcare for replacement through the life of your Corometrics monitor. The combo transducers, for both ultrasound and toco, feature the same trusted technology—but now with a novel combination of button and loop styles in a single transducer.

There is no substitute for GE Healthcare combo transducers. Here's why:



Reliability

- Verified compatible with Corometrics 250cx series and 170 series monitors by GE Healthcare Corometrics engineers
- Proven Corometrics technology with a 50-year track record



Simplicity

- Convenient support thanks to a single-source solution for monitors and transducers
- Combo button-and-loop design provides flexibility in a single transducer



Parameter excellence

- Our nine-crystal technology helps achieve a wide focal region for more uniform coverage at greater depths than conventional seven-crystal transducers
- Designed to support minimal repositioning even with challenging patients of varying sizes, shapes, and gestational ages



Patient comfort

- Flat toco transducer designed for increased comfort
- Watertight for cleanability and in-water use



There are three ways you can order combo transducers:

- **Phone:** GE Customer Service
1-800-437-1171, Option 2, Option 3
- **Online:** GE Service Shop
services.gehealthcare.com/gehcstorefront
- **Email:** GEHCParts.OrderInquiry@ge.com

Part number	Description
2108346-001	Ultrasound transducer combo
2108347-001	Toco transducer combo



[gehealthcare.com](https://www.gehealthcare.com)

© 2019 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Corometrics are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB66326US