



Fetal Spiral Electrodes and Cables

Quality accessories. Accurate measurements.

In labor and delivery, patient comfort and safety are top priorities. Our fetal spiral electrodes and related products are designed to help you deliver high-quality care by providing accurate measurements.

The electrodes and cables are designed to be easy and convenient to use. The electrodes have lead wires fully protected by a safety cap. A self-adhesive grounding electrode attaches to the patient.

Features

- Drive and guide tube removable in one step
- Guide tube flexible to enable gentle shaping for ease of insertion
- Simplified leg plate connection requiring no orientation
- Spiral recessed to guard against needle sticks when opening the package

Product compatibility

The accessories listed below are compatible and verified with the following monitoring devices:

- Corometrics™ 174 fetal monitor
- Corometrics 259cx fetal monitor
- GE Mini Telemetry System

Our equipment. Our accessories. Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step

PRODUCT SPECIFICATIONS

Materials

Not made with natural rubber latex or DEHP
ROHS (Restriction of hazardous substances) compliant

Operating conditions

Recommended ambient temperature 10° to 35°C (50° to 95°F)

Storage conditions

Recommended temperature -10° to 50°C (14° to 122°F)

Cleaning solutions (cable only)

Ethanol alcohol maximum	96% by volume
Isopropyl alcohol maximum	60% by weight
Glutaraldehyde	2%
Phenol	2%

The electrode and the pad are single use and the cable is re-usable.

The cable is not intended to be sterilized or autoclaved.

Compliance

CE Mark

PART LIST

Part Number	Description	Quantity
2109004-001	31479812, Fetal Spiral Electrode Cable, Round connector	1
2109002-001	31479549, Fetal Spiral Electrode, Single Helix, Single Patient Use	50
2109006-001	50000095, Fetal Spiral Electrode, Attachment Pad, Large, Single Patient Use	50
2109005-001	50000107, Fetal Spiral Electrode Cable, Rectangular connector	1

Product may not be available in all countries and regions.

Contact a GE Healthcare Representative for more information.

Please visit www.gehealthcare.com/promotional-locations



Manufactured by:

Cardinal Health 200, LLC
7000 Cardinal Place
Dublin, OH 43017
USA



© 2019 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Corometrics are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2146773 1/19