




# GE ECG Electrodes



## Secure adhesion. Quality signal transmission.

GE electrodes for patient monitoring and telemetry are designed to ensure positive contact with the skin and enable reliable connections with ECG leadwires. A broad family of electrodes includes configurations to accommodate a full range of patients and a variety of skin conditions.

- Variety of styles, including choice of snaps or tabs
- Adult, pediatric and neonate sizes available
- Designed for easy placement and removal
- Choice of cloth, foam or tape; solid or wet gel
- Round, square, rectangle and teardrop shapes
- Hypoallergenic materials; gentle adhesives available for sensitive skin
- Low profile to prevent accidental detachment by bed sheet or gown
- Resealable pouches with special lining to prevent electrodes from drying out
- Water-based electrode gel to support high signal integrity

## Product compatibility

GE ECG electrodes are compatible with GE ECG monitoring, telemetry and diagnostic resting equipment.

## Accessories simply compatible with your equipment.

You can trust GE Healthcare accessories and supplies to keep your clinical devices working efficiently.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step

## PRODUCT SPECIFICATIONS

### Materials

---

Not made with natural rubber latex, DEHP or BPA  
RoHS (Restriction of hazardous substances) compliant

### Storage conditions

---

Recommended temperature 10° to 32°C (50° to 90°F).  
Electrodes should be stored in the original pouch.

### Shelf life data

---

Varies between 18 to 24 months from the manufacturing date. See details on back.

### Certifications

---

Complies with AAMI Disposable ECG electrode standards

## Part List

Part number	Description	Population	Conductive Medium	Material				Form	Dimensions	Quantity/ Pouch	Quantity/ Case
				Foam	Cloth	Tape	Adhesive				
2106924-001	ECG Snap* Electrode	Adult	Conductive Solid Hydrogel	X				Teardrop	38.1mm/1.5in	30	600
2106890-001	ECG Snap* Electrode	Adult	Conductive Solid Hydrogel		X			Teardrop	47.625mm/1.875in	60	600
2106889-001	ECG Snap* Electrode Radiotranslucent	Adult	Conductive Solid Hydrogel	X				Teardrop	47.625mm/1.875in	60	600
2106915-001	ECG Snap* Electrode	Adult	Conductive Solid Hydrogel	X				Teardrop	47.625mm/1.875in	60	600
2104779-001	ECG Snap** Electrode	Adult	Conducive 10% KCl Wet gel	X				Teardrop	38.1mm/1.5in	60	600
2104777-001	ECG Snap* Electrode	Adult	Conductive Solid Hydrogel			X		Teardrop	47.625mm/1.875in	60	600
2106921-001	ECG Snap** Electrode	Adult	Conductive 7% KCl Wet Gel	X				Teardrop	38.1mm/1.5in	5/strip 60/Pouch	600
2104780-001	ECG Tab* Electrode	Adult	Conductive Solid Hydrogel				X	Large Tab	28.575mm x 36.525mm/ 1.125 x 1.438in	10/Strip 100/Pouch	4000 (4 boxes of 1000)
2106911-001	ECG Tab* Electrode	Adult	Conductive Solid Hydrogel				X	Medium Tab	34.036mm x 22.86mm/ 1.34in x .90in	10/Strip 100/Pouch	4000 (4 boxes of 1000)
2106917-001	ECG Snap* Electrode	Pediatric	Conductive Solid Hydrogel		X			Rectangle	22.225mm x 33.34mm/ 0.875in x 1.3125in	3/Strip 60/Pouch	600
2106918-001	ECG Snap* Electrode	Pediatric	Conductive Solid Hydrogel	X				Round	25.4mm/1in	3/Strip 60/Pouch	600
2106919-001	ECG Snap* Electrode Radiotranslucent	Pediatric	Conductive Adhesive Hydrogel		X			Square	22.225mm/0.875in	3/Strip 60/Pouch	600
2104778-001	ECG Tab* Electrode	Pediatric	Conductive Adhesive Hydrogel				X	Small Tab	12.7mm x 32.54mm/ 0.5in x 1.281in	10/Strip 100/Pouch	4000 (4 boxes of 1000)
2106891-001	ECG Pre-wired* AHA Electrode	Neonatal	Conductive Adhesive Hydrogel		X			Square	22.225mm/0.875in	3	300 (100 sets)
2106892-001	ECG Pre-wired* IEC Electrode	Neonatal	Conductive Adhesive Hydrogel		X			Square	22.225mm/0.875in	3	300 (100 sets)

\* = shelf life 24 months | \*\* = shelf life 18 months

**Product may not be available in all countries and regions.**  
Contact a GE Healthcare Representative for more information.  
[www.gehealthcare.com/promotional-locations](http://www.gehealthcare.com/promotional-locations)

Distributed by:  
GE Healthcare Finland Oy



Manufactured by:  
Graphic Controls  
400 Exchange Street  
Buffalo, NY 14204  
USA

© 2018 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram, are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2146750 rev2 11/18