



AHA compliant

Reusable ECG Leadwires

Easy connection. Accurate measurements.

GE ECG leadwires help you reliably make the correct connection to GE monitoring devices. These high-quality reusable leadwires complement your monitoring equipment, enabling accurate measurements.

GE leadwires undergo strict testing to help ensure quality, reliability and safety. They have been tested and verified as compatible with GE monitoring devices¹ to deliver optimum performance and a high standard of care.



Durable

Chemical-resistant jacketing extends service life, facilitates cleaning and remains flexible to withstand bending and twisting and enhance patient comfort.

Ease of use and efficient transfer

The leadwires are fast and easy to apply with clearly labeled, color-coded snaps and grabbers. Reusable leadwires can help improve efficiency in patient transfers by allowing patients to wear only one dedicated leadwire set throughout their hospital stay, if the facility is completely standardized.

Features

GE ECG lead wires are available in 3-, 5-, 6- and V-lead configurations. Leadwires are electrically shielded to protect ECG signal integrity.

Product compatibility

The ECG leadwires are compatible and validated with the following GE equipment:

- CARESCAPE™ Monitors B450/B650/B850
- CARESCAPE Patient Data Module (PDM) and Patient Side Module (PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Solar™ monitors with CARESCAPE Patient Data Module
- ApexPro™ CH telemetry

ECG trunk cables are also compatible with GE ECG leadwire sets.

Accessories simply compatible with your equipment.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step

PRODUCT SPECIFICATIONS

Materials

Does not contain natural rubber latex, or DEHP
ROHS (Restriction of hazardous substances) compliant

Operating conditions

Recommended
ambient temperature 10° to 35°C (50° to 95°F)

Storage conditions

Recommended
temperature -10° to 50°C (14° to 122°F)

Cleaning solutions

Ethyl alcohol maximum	90% by volume
Isopropyl alcohol maximum	60% by weight
Glutaraldehyde	2%
Phenol	2%

Do not sterilize or autoclave.

Compliance

American Heart Association (AHA) standards
CE Mark



ORDERING INFORMATION

Part Number	Description
3-Lead (Grabber)	
2106390-001	ECG Leadwire set, 3-lead grouped, Grabber, AHA, 74 cm/29 in
2106390-002	ECG Leadwire set, 3-lead grouped, Grabber, AHA, 10 cm/ 51 in
5-Lead	
2106391-001	ECG Leadwire set, 5 -lead, Grabber, AHA, 74 cm/ 29 in
2106391-002	ECG Leadwire set, 5-lead, Grabber, AHA, 130 cm/ 51 in
2106391-005	ECG Leadwire set, 5-lead, Grabber, AHA, mix 74 cm/ 29 in, 130 cm/51 in
2106393-001	ECG Leadwire set, 5-lead V2-6, Grabber, AHA, 130 cm/51 in
2106393-003	ECG Leadwire set, 5-lead V2-6, Grabber, AHA, 74 cm/29 in
2106389-001	ECG Leadwire set, 5-lead grouped, Grabber, AHA, 74 cm/ 29 in
2106389-002	ECG Leadwire set, 5-lead grouped, Grabber, AHA, 130 cm/51 in
2106389-005	ECG Leadwire set, 5-lead grouped, Grabber, AHA, mix 74 cm/ 29 in, 130 cm/ 51 in
6-Lead	
2106397-001	ECG Leadwire set, 6-lead, Grabber, AHA, 74 cm/ 29 in
2106397-002	ECG Leadwire set, 6-lead, Grabber, AHA, 130 cm/ 51 in
2106397-003	ECG Leadwire set, 6-lead, Grabber, AHA, mix 74 cm/ 29 in, 130 cm/51 in

Part Number	Description
3-Lead (Snap)	
2106385-001	ECG Leadwire set, 3-lead, snap, AHA, 74 cm/ 29 in
2106385 -002	ECG Leadwire set, 3-lead, snap, AHA, 130 cm/ 51 in
5-Lead	
2106381-001	ECG Leadwire set, 5-lead, snap, AHA, 74 cm/ 29 in
2106381-002	ECG Leadwire set, 5-lead, snap, AHA, 130 cm/ 51 in
2106381-005	ECG Leadwire set, 5-lead, snap, AHA, mix 74 cm/ 29 in, 130 cm/ 51 in
2106383-001	ECG Leadwire set, 5-lead, grouped, snap, AHA, 74 cm/ 29 in
2106383-002	ECG Leadwire set, 5-lead, grouped, snap, AHA, 130 cm/ 51 in
2106383-005	ECG Leadwire set, 5-lead, grouped, snap, AHA, mix 74 cm/ 29 in, 130 cm/ 51 in
2106394-001	ECG Leadwire set, 5-lead V2-6, snap, AHA, 130 cm/ 51 in
2106394-003	ECG Leadwire set, 5-lead V2-6, snap, AHA, 74 cm/ 29 in
6-Lead	
2106395-001	ECG Leadwire set, 6-lead, grouped, snap, AHA, 74 cm/ 29 in
2106395-002	ECG Leadwire set, 6-lead, grouped, snap, AHA, 130 cm/ 51 in
2106395-003	ECG Leadwire set, 6-lead, grouped, snap, AHA, mix 74 cm/ 29 in, 130 cm/ 51 in
Leadwire Separator	
2106399-001	Leadwire Separator, 3 pieces

Product may not be available in all countries and regions.
Contact a GE Healthcare Representative for more information.
Please visit www.gehealthcare.com/promotional-locations

Imagination at work

© 2018 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram, are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE, the GE monogram, Imagination at work, CARESCAPE, ApexPro CH and Solar are trademarks of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2145248 6/18