

Comparing GE HealthCare's Resting ECG Systems

Get more out of your resting ECG with the MAC™ family. Discover which ECG is right for your facility.



	MAC 5	MAC 7*	MAC VU360™
<p>● Standard feature — Not available ▲ Optional feature</p>			
Display capability			
Color	●	●	●
12 waveforms	●	●	●
15 waveforms	—	▲	●
Waveform quality check (hook-up advisor)	●	●	●
Auto ECG upon clean hookup	▲	▲	●
Touchscreen	●	●	●
Size	8.9" 892 x 558	10.1" 1280 x 800	15.6" Full HD
Analysis programs			
Measurements	●	●	●
Interpretation with age and gender software	▲	●	●
Left Ventricular Hypertrophy	●	●	●
High definition pacemaker detection	●	●	●
Right Ventricular Involvement (RVI) 15-lead	—	▲	●

● Standard feature — Not available △ Optional feature

	MAC 5	MAC 7	MAC VU360
Analysis programs cont.			
Vectorcardiography	—	△	●
Acute Coronary Syndrome (ACS)	△	△	△
Critical values — user definable	△	△	△
Prior ECGs	—	△	△
Digital rhythm	△	△	●
Full disclosure	△ 5 mins	△ 5 mins	△ 5 mins
Communication			
Wireless LAN	△	△	●
Wired LAN	●	●	●
Memory/Data storage			
Internal storage	300 records consisting of 10 second Resting ECG and rhythm records; 200 minutes Digital Rhythm or Full Disclosure records	1000 records consisting of 10 second Resting ECG records and Digital Rhythm records on the device internal memory	1000 records consisting of 10 second Resting ECG records and Digital Rhythm records on the device internal memory
External memory — SD card / USB flash	●	●	●
Training			
Application Training Kits	—	△	△
Workflow			
MUSE™** ECG Management System	●	●	●
CardioSoft™ ECG Management System (transmit ECG records)	●	●	●
Orders	△	△	●
ADT	△	△	●
Integrated EMR Connectivity ¹	—	—	△
DICOM ^{®2} /EMR Gateway Connectivity	△	△	△
Shared folder / FTPS export	●	●	●
Barcode	△	△	△
PDF export	●	●	●
XML export	△	△	△
Web enabled	—	—	△
Network printing	△	△	△
ReadyFix™ remote fleet management ³	—	—	△
Paper size			
Thermal Z-fold paper	A4/Letter or A5	A4/Letter	A4/Letter
Single package size	150 Sheets	150 Sheets	150 or 300 Sheets
Design			
Height adjustable trolley	—	—	△
Enclosure storage	—	—	△
Designed and tested for cleaning/disinfection	●	●	●

● Standard feature — Not available ▲ Optional feature

Hot swappable batteries, extra battery charge options

	MAC 5	MAC 7	MAC VU360
	—	—	●
Specifications	A4/Letter or A5 or No printer		
Height	12.5 / 12 / 9.5 cm	22 cm	145 cm
Width	37 / 27 / 27 cm	33 cm	52 cm
Depth	32.5 / 32.5 / 32 cm	40.5 cm	60 cm
Weight (including battery, without paper)	4 / 3.4 / 2.3 kg	4.8 kg	34 kg

Augmented reality

Our augmented reality experience allows you to explore key features and see how easily MAC ECG products can fit into your existing office or hospital environment

Point your camera at the QR code. Tap the banner that appears on your screen.



MAC 5 A4



MAC 5 A5



MAC 5 Lite



MAC 7



MAC VU360

Please contact your GE HealthCare representative to learn more, or visit www.gehealthcare.com

* MAC 7 version 2.

** Customer provided display (supported range).

1. Via MUSE (v8 or higher).
2. Via MUSE (v9 or higher).
3. Additional subscription costs apply. MAC VU360 devices must be updated to v3 software to enable ReadyFix.

© 2025 GE HealthCare.

Not all products or features are available in all markets. Contact a GE HealthCare Representative for more information. Please visit <http://www.gehealthcare.com>. Data subject to change. MAC, MAC VU360, CardioSoft, ReadyFix, and MUSE are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license. DICOM is a registered trademark of National Electrical Manufacturers Association (NEMA).

JB03331XU | December 2025



GE HealthCare