

Optima™ AMX Navigate, Fixed Column

GoldSeal™ Refurbished

The GoldSeal AMX Navigate is a revolutionary mobile X-ray designed to have you in control, position effortlessly and take and process images flawlessly. Make key decisions smart and accelerate your course of action. It delivers consistent, reliable results whenever and wherever you need it.

Features

- Helix Advanced Image Processing 2.1: Delivers extraordinary anatomical detail and AI-based brightness and contrast adjustment
- FlashPad™ HD Detectors: 100-micron pixel pitch, 75% DQE offering high image clarity. Multiple sizes available
- QuickShare: FlashPad HD detectors can be shared across compatible GE HealthCare Fixed or Mobile X-Ray systems
- QuickConnect: Adaptive wireless technology to improve image transfer
- QuickEnhance: Simple 1-touch image reprocessing to support anatomy specific customized visualization
- QuickClean: The screen itself is easy to sanitize due to the QuickClean touchscreen lockout application
- Touchscreen monitor: 21.5 in touchscreen with intuitive user interface
- Generator: 30 kW
- Column: The system offers a fixed column (does not include the power-assisted telescoping column)
- IT Security Package: Delivers cybersecurity features which can be configured as needed



Available options

- Zero Click Exam: Helps streamline the imaging process. Zero Click Exam leverages RFID Badge Swipe, barcode patient verification, automated protocol selection and AI processing
- Accessories: Grids, Detector Grips, Battery options
- Optional software: AutoGrid™ and Critical Care Suite
- HIS/RIS integration
- Onsite clinical applications training and TiP™ Virtual Assist

Warranty

Includes full 1 year warranty on the unit. detector coverage information is identified in the GE HealthCare Warranty statement.

NEMA XR 29-2013 standard

This product is compliant with the NEMA XR 29-2013 standard.

Contact us

GoldSeal systems are quoted subject to availability. To confirm current availability of systems or additional options in your region, contact your GE HealthCare representative.

Your access to quality medical imaging is our priority. Get the most for your budget by choosing GoldSeal today. If it meets our standards, we know it will meet yours.

Learn more at www.gehealthcare.com/goldseal

Good Refurbishment Practice

The GoldSeal program is consistent with Good Refurbishment Practice (GRP) established by leading industry trade associations such as MITA (Medical Imaging & Technology Alliance) a division of NEMA (National Electrical Manufacturers Association), COCIR (European Coordination Committee of the Radiological, Electromedical and Healthcare IT Industry) and JIRA (Japan Industries Association of Radiological Systems).