



GE HealthCare

LOGIQ e X[◇]

Premium, AI-powered console-quality ultrasound in a compact system for anywhere care is delivered*



For Rheumatologists

LOGIQ™ e X brings the imaging performance of a console system into a lightweight, laptop-sized ultrasound built for rheumatologists to deliver consistent, objective care. With advanced imaging capabilities, AI-powered tools and workflow-enhancing features, LOGIQ e X helps rheumatologists track disease progression and treat toward remission with confidence and reproducibility. LOGIQ e X empowers clinicians with simplicity, speed and confidence to deliver high-quality care in today's most dynamic environments.

Simple

Easy to learn and operate



Cart options: Two compact cart options for comfortable, mobile scanning: one with height adjustment and extended battery and one with a fixed-height design for streamlined simplicity



Enhanced ergonomics: Features laptop-style portability, dual battery and intuitive control layout for comfortable scanning

Fast

Designed to keep pace with busy environments



Voice Control: Activates system functions using simple voice commands† to enhance workflows, reduce keystrokes and support ergonomics



Scan Assistant: Customized scanning program that guides users step by step through exam protocols which can reduce keystrokes and shorten exam time

Confident

Advanced visualization tools and AI-powered features that help support confident care




Speckle Reduction Imaging (SRI) – HD: Enhances visibility of organs and lesions with high-definition contrast resolution that suppresses speckle artifacts while preserving true tissue architecture

New



MSK Go

Highlights musculoskeletal anatomical structures on acquired ultrasound images to support clinician confidence and help streamline workflow and reduce exam time



More capable. More versatile. More portable.

The LOGIQ e X delivers simplicity, speed and confidence in a compact laptop design.

Large 15.6" screen

New color

Trackball/trackpad options


Intuitive control layout

Dual battery

(Up to 90 min. scan time on console;
up to 3.5 hours scan time on cart)

Two-/three-probe
port connector

Cart options



Ideal for:

Rheumatology

Radiology

Orthopedics

Rehabilitation

Physiotherapy

Compatible probes



Linear array:

L10-22-RS

L8-18i-RS

ML6-15-RS

L4-12t-RS

9L-RS

L6-12-RS

12L-RS

L4-20t-RS

Convex array:

C1-5-RS

C2-9-RS

4C-RS

Microconvex array:

8C-RS

Wireless:


Vscan Air™ CL (curved/linear)

Vscan Air SL (sector/linear)


Powered by XDclear™ technology

Clear images. Consistent exams. Confident results.


Advanced imaging capabilities deliver exceptional clarity to make informed decisions.



MSK Go



SRI-HD



Power Doppler Imaging (PDI)

◇ LOGIQ e X is a marketing representation of the LOGIQ e product; feature availability may vary depending on region.

* The device is intended for use in professional healthcare facilities. Use is restricted to environmental properties described in the user manual.

† Voice Control currently supports English, Chinese, Latin American Spanish, Swedish, German, Italian, French, Russian, Greek, Brazilian Portuguese, Dutch and Japanese.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with your local GE HealthCare representative for details.

© 2026 GE HealthCare. LOGIQ and XDclear are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

June 2026
JB37469XX



GE HealthCare