



GE HealthCare

LOGIQ e X[◇]

Premium console-quality ultrasound in a compact system for anywhere care is delivered*



For Vascular care

LOGIQ™ e X brings the imaging performance of a console system into a lightweight, laptop-sized ultrasound built for providing precise vascular care. With advanced imaging capabilities and workflow-enhancing features, LOGIQ e X provides detailed vascular visualization and needle guidance to support accurate placement during a variety of procedures. LOGIQ e X empowers clinicians with simplicity, speed and confidence to deliver high-quality care in today's most dynamic environments.

Simple

Easy to learn and operate



Cart options: Two compact cart options for comfortable, mobile scanning: one with height adjustment and extended battery and one with a fixed-height design for streamlined simplicity



Enhanced ergonomics: Features laptop-style portability, dual battery and intuitive control layout for comfortable scanning

Fast

Designed to keep pace with busy environments



Voice Control: Activates system functions using simple voice commands† to enhance workflows, reduce keystrokes and support ergonomics



Scan Assistant: Customized scanning program that guides users step by step through exam protocols which can reduce keystrokes and shorten exam time

Confident

Advanced visualization tools that help support confident care



Needle Recognition: Designed to enhance the visibility and accuracy of needle placement during ultrasound-guided procedures by enhancing needle visualization in real-time

New



B-Flow

Delivers real-time, high-resolution blood flow visualization and provides clear, angle-independent visualization of small vessels for confident assessment of flow dynamics and vascular lesions



More capable. More versatile. More portable.

The LOGIQ e X delivers simplicity, speed and confidence in a compact laptop design.

Large 15.6" screen

New color

Trackball/trackpad options


Intuitive control layout

Dual battery

(Up to 90 min. scan time on console;
up to 3.5 hours scan time on cart)

Two-/three-probe
port connector

Cart options



Ideal for:

Vascular clinics

Surgery centers

Vein clinics

Dialysis clinics

Compatible probes



Linear array:

L4-20t-RS

L8-18i-RS

ML6-15-RS

L4-12t-RS

9L-RS

L6-12-RS

12L-RS

L10-22-RS

Phased array:

3Sc-RS

6S-RS

Convex array:

C2-9-RS

C1-5-RS

4C-RS

Microconvex array:

8C-RS

Wireless:


Vscan Air™ CL (curved/linear)

Vscan Air SL (sector/linear)


✦ Powered by XDclear™ technology

Clear images. Consistent exams. Confident results.

Advanced imaging capabilities deliver exceptional clarity to make informed decisions.



Scan Assistant



B-Flow



Needle Recognition

✦ LOGIQ e X is a marketing representation of the LOGIQ e product; feature availability may vary depending on region.

* The device is intended for use in professional healthcare facilities. Use is restricted to environmental properties described in the user manual.

† Voice Control currently supports English, Chinese, Latin American Spanish, Swedish, German, Italian, French, Russian, Greek, Brazilian Portuguese, Dutch and Japanese.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with your local GE Healthcare representative for details.

© 2026 GE Healthcare. LOGIQ, Vscan Air, and XDclear are trademarks of GE Healthcare. GE is a trademark of General Electric Company used under trademark license.

May 2026
JB37475XX



GE Healthcare