



GE HealthCare

LOGIQ e X[◇]

Premium, AI-powered console-quality ultrasound in a compact system for anywhere care is delivered*



For Veterinary teams

LOGIQ™ e X brings the imaging performance of a console system into a lightweight, laptop-sized ultrasound built for veterinary clinics and teams that travel between clinics. With advanced imaging capabilities, workflow-enhancing features and broad probe compatibility, LOGIQ e X helps provide diagnostic confidence across different animal types. LOGIQ e X empowers veterinary teams with simplicity, speed and confidence to deliver high-quality care in today's most dynamic environments.

Simple

Easy to learn and operate



Two-three probe port connector: Enables the connection of two or three ultrasound probes to a single system allowing seamless switching between probes



Enhanced ergonomics: Features laptop-style portability, dual battery and intuitive control layout for comfortable scanning

Fast

Designed to keep pace with busy environments



Voice Control: Activates system functions using simple voice commands† to enhance workflows, reduce keystrokes and support ergonomics



Scan Assistant: Customized scanning program that guides users step by step through exam protocols which can reduce keystrokes and shorten exam time

Confident

Advanced visualization tools that help support confident care



Speckle Reduction Imaging (SRI) – HD: Enhances visibility of organs and lesions with high-definition contrast resolution that suppresses speckle artifacts while preserving true tissue architecture



More capable. More versatile. More portable.

The LOGIQ e X delivers simplicity, speed and confidence in a compact laptop design.

Large 15.6" screen

New color

Trackball/trackpad options


Intuitive control layout

Dual battery

(Up to 90 min. scan time on console;
up to 3.5 hours scan time on cart)

Two-/three-probe
port connector

Cart options



Ideal for:

Equine

Companion animals

Farm animals

Wild and exotic
animals

Compatible probes



Linear array:

L10-22-RS

L8-18i-RS

ML6-15-RS

L4-12t-RS

9L-RS

L6-12-RS

12L-RS

L4-20t-RS

LK760-RS

Phased array:

3Sc-RS

12S-RS

6Tc-RS

6S-RS

Convex array:

C1-5-RS

C2-9-RS

4C-RS


Microconvex array:

8C-RS

E8C-RS

IC9-RS

✦ Powered by XDclear™ technology



SRI-HD

Clear images. Consistent exams. Confident results.

Advanced imaging capabilities deliver exceptional clarity to make informed decisions.

◇ LOGIQ e X is a marketing representation of the LOGIQ e product; feature availability may vary depending on region.

* The device is intended for use in professional healthcare facilities. Use is restricted to environmental properties described in the user manual.

† Voice Control currently supports English, Chinese, Latin American Spanish, Swedish, German, Italian, French, Russian, Greek, Brazilian Portuguese, Dutch and Japanese.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with your local GE HealthCare representative for details.

© 2026 GE HealthCare. LOGIQ, and XDclear are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

May 2026
JB37474XX



GE HealthCare