



# CRITIKON SOFT-CUF™

## Blood Pressure Cuffs



Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer-term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It can be cleaned and reused on multiple patients, but is priced economically so it can be disposed of after a single patient use, helping you affordably prevent infections.

### Designed for patient safety, comfort and accuracy

- Limited reuse, yet economical enough for pairing with a single patient throughout their hospital stay to help with your infection control efforts
- Soft, absorbent material is not made with DEHP or natural rubber latex to help enhance patient safety
- Designed with rounded corners for patient comfort
- Designed and clinically validated for performance accuracy to help you make well-informed decisions
- DINACLICK™ connector meets international standards and helps prevent misconnections to luer fittings
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit

### Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors and a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility









### Cuff assortment packs available

- For low-volume cuff utilizers, three types of assortment packs are available for easier ordering and more optimal inventory management



DINACLICK  
Connector



SOFT-CUF			GE Healthcare				Non-GE Healthcare			
										
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Bayonet Connector	1-Tube Submin Connector	SOFT-CUF Inflation System
Infant	Orange	8-13 cm	2500	SFT-P1-2A	2401	2450 <sup>1</sup>	002331	2360	2320	
Child	Green	12-19 cm	2501	SFT-P2-2A <sup>2</sup>	2402	2451	002332	2361	2321	002584
Child Long	Green	12-19 cm	2506	SFT-P2-2A-L	2400					
Small Adult	Lt. Blue	17-25 cm	2502 <sup>2</sup>	SFT-A1-2A <sup>2</sup>	2403 <sup>2</sup>	2452	002333	2362 <sup>2</sup>	2322	002583
Small Adult Long	Lt. Blue	17-25 cm	2507	SFT-A1-2A-L	2407		002334			
Adult	Navy	23-33 cm	2503 <sup>2</sup>	SFT-A2-2A <sup>2</sup>	2404 <sup>2</sup>	2453	002335	2363 <sup>2</sup>	2323	002580
Adult Long	Navy	23-33 cm	2604 <sup>2</sup>	SFT-A2-2A-L <sup>2</sup>	2116 <sup>2</sup>	2454	002336	2140 <sup>2</sup>		002582
Large Adult	Rose	31-40 cm	2504 <sup>2</sup>	SFT-A3-2A <sup>2</sup>	2405 <sup>2</sup>	2455	002337	2364 <sup>2</sup>	2324	002581
Large Adult Long	Rose	31-40 cm	2644	SFT-A3-2A-L	2117	2456	002338	2141		002585
Thigh	Brown	38-50 cm	2505	SFT-T1-2A	2406	2457	002339	2365	2325	

1. Infant BP Cuff for Datex-Ohmeda\* monitors has male slip luer connector

2. Has a respective part number with an "-INT" suffix, denoting that the Instructions For Use (IFU) is available in multiple languages

Note: 20 cuffs/box

SOFT-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	002695	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
2-Tube Submin Connectors	2056304-002	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2298	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
2-Tube DINACLICK Connectors	2059303-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)



### Imagination at work

[www.gehealthcare.com](http://www.gehealthcare.com). Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Data subject to change.

© 2016 General Electric Company.

GE, GE Monogram, Imagination at work, CRITIKON, SOFT-CUF, and DINALICK are trademarks of General Electric Company.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.