



GE Healthcare

Clinical Accessories that support **Infection Prevention & Control**

Infection Prevention and Control Policy is to ensure that the Hospital's provide an environment and system of care which minimises the risk of infection to patients, staff and visitors.

GE Accessories offer a reliable disposable solution that delivers value.

- Single patient use accessories **reduce the risk of cross-contamination** and help **lower the risk of costly equipment associated infections**
- No sterilization or assembly required, which could lead to **time and cost savings**
- **Help reduce the risk of nosocomial infections**

INTENSIVE CARE

Single-Patient-Use Valve



Flow Sensor Only
1505-3231-001



Housing only
(includes water trap)
1505-8568-001



Whole assembly
(Housing and Flow sensor)
1505-3848-001



Inspiratory Safety Guard
2083208-001 (10 pcs/Pkg)

NIBP cuffs



Yellow Critikon™ Classic-Cuf*



Soft-Cuf*



Radial-Cuf*

Disposable SpO₂ Sensor



TruSignal AllFit Sensor
TS-AF-10 (10 pcs/Pkg)
TS-AF-25 (25 pcs/Pkg)



TruSignal Adult/Pediatric Sensor
TS-AP-10 (10 pcs/Pkg)
TS-AP-25 (25 pcs/Pkg)



D-fend+
10/case
881319-HEL¹



D-fend Pro+
10/case
M1200227¹



Gas Sampling line
2m - 10/case
2097307-001

Airway Gas Sampling Set, Disposable Adult/Pediatric, TPE, 2m/10ft
20 pcs/Pkg
8090121521

D-lite+ sensor
50/case
896952

PERIOPERATIVE CARE

Gas Sampling set



Gas Sampling line
2m - 10/case
2097307-001

Airway Gas Sampling Set, Disposable Adult/Pediatric, TPE, 2m/10ft
20 pcs/Pkg
8090121521

D-lite+ sensor
50/case
896952

Water trap



D-fend
10/case
876446-HEL²



D-fend Pro
10/case
M1182629²

HEPA Filter



Adult HEPA pleated paper hydrophobic bacterial viral filter and Luer port
2106570-008

NIBP cuffs

SOFT-CUF			GE Healthcare				Non-GE
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Bayonet Connector
Small Adult	Light Blue	17-25 cm	2502 ³	SFT-A1-2A ³	2403 ³	2452	2362 ³
Small Adult Long	Light Blue	17-25 cm	2507	SFT-A1-2A-L	2407		
Adult	Navy	23-33 cm	2503 ³	SFT-A2-2A ³	2404 ³	2453	2363 ³
Adult Long	Navy	23-33 cm	2604 ³	SFT-A2-2A-L ³	2116 ³	2454	2140 ³
Large Adult	Rose	31-40 cm	2504 ³	SFT-A3-2A ³	2405 ³	2455	2364 ³

RADIAL-CUF		GE Healthcare			
Color	Forearm Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector
Purple	26-36 cm ⁴	SFT-F1-2E	SFT-F1-2A	SFT-F1-2C	SFT-F1-2D

Note: GE Critikon RADIAL-CUF BP cuffs come 20 in a box.

YELLOW CLASSIC-CUF			GE Healthcare		
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector
Small Adult	Yellow/Blue	17-25 cm	2607	CLA-A1-2A-X	2672
Adult	Yellow/Blue	23-33 cm	2602	CLA-A2-2A-X	2673
Large Adult	Yellow/Blue	31-40 cm	2642	CLA-A3-2A-X	2674

GE masks and breathing circuits are also available at your service, please contact your local sales representatives for further information.

gehealthcare.com

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information.

GE, GE Monogram are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company.

¹ Recommended to change water trap for each new patient, whenever the host device prompts you to do so, or every 24 hours.

² Change whenever the host device prompts you to do so. The maximum lifetime of a water trap is two months.

³ Has a respective part number with an "-INT" suffix, denoting that the Instructions for Use (IFU) is available in multiple languages.

⁴ Forearm measurement