



Mac-Lab/CardioLab v6.9.6 Virtual Review Software Support

Product Group: Interventional Invasive Products

Products: Mac-Lab IT/XT/XTi, CardioLab IT/XT/XTi,
SpecialsLab and ComboLab IT/XT/XTi
Recording Systems

Subject: Virtual Review Supported Software Versions

Date: Sept 8, 2015

Summary

This document details the latest software compatibility for the Virtual Review. Where there are differences between the product manual and this document, this document supersedes information provided in the manual.

v6.9.6 Virtual Review Software Support

- **VMWare Player:** VMware Player 7.0.0 Pro or VMware Player 7.1.0 Pro (recommended). Other versions have not been tested for compatibility.

NOTE: The Virtual Review is not compatible with the VMware Player's Unity mode.

NOTE: VMware has revised and renamed Player Pro 7.x product to Workstation Player 12. VMware Player Pro 7.0.0/7.1.0 license and binaries are available by purchasing Workstation Player 12 and downgrading the license within the VMware's MY VMware functionality on the VMware website. The following URL provides instruction on how to downgrade any VMware license:

http://kb.vmware.com/selfservice/microsites/search.do?language=en_US&cmd=displayKC&externalId=2006975