


How the best see better

Introducing TrueFidelity™ CT Images from GE Healthcare





Vision that pushes radiology further.

Artificial Intelligence. Machine Learning. Deep Learning. Radiology is undergoing a transformation.

Deep Learning Image Reconstruction (DLIR) promises unparalleled benefits for patients, along with the radiologists and technologists dedicated to their care. And after nearly half a century at the forefront of computed tomography, GE Healthcare is uniquely positioned to ensure this latest advance keeps its promise.


GE Healthcare pioneered and consistently pushed the science of image reconstruction further. TrueFidelity CT Images are more than a radical, next generation improvement. They elevate the vision of what you and Deep Learning Image Reconstruction can achieve—together.


Where Deep Learning does its learning matters.

A Deep Learning Image Reconstruction application is only as good as the training it receives.

GE Healthcare's proprietary DLIR training reflects our unmatched understanding of what successful DLIR requires. And what radiologists want. The foundation of that training is GE's library of thousands of low noise, filtered back projection images.

These ground truth images cover every anatomy and are the gold standard for image quality. The artificial intelligence that powers our Deep Learning Image Reconstruction gets its education here. An education that allows you and TrueFidelity CT Images to achieve never-before-possible clarity at low dose.





In supervised training, GE's Deep Neural Network (DNN)-based DLIR engine compares the output to thousands of ground truth images.


Confidence. Not compromise.

Compared with even the most sophisticated Model-Based Iterative Reconstruction, TrueFidelity CT Images are scanning taken to another level.


Contrast visualization is maintained; noise and artifacts are minimized; edges are maintained—just enough—so there's remarkable clarity and none of the compromise that comes with unfamiliar noise texture.*

The result is an easy, information-rich interpretation experience. An experience that gives diagnosticians the confidence they require, even as it potentially improves scan read times and fights radiologist fatigue.


*DLIR neural networks give an image appearance (as shown on NPS plots) similar to traditional high-dose, low-noise FBP images. Demonstrated in phantom testing comparing images reconstructed from the same raw data, with DLIR-H and ASiR-V 100%, using the standard kernel.






0.625 mm FBP



0.625 mm ASiR-V 50%



0.625 mm TrueFidelity



TrueFidelity

“This is what a head CT image should look like and never did, even in the high-dose era. Less noise, fewer artifacts, and a pleasure to interpret.”

—**Lawrence N. Tanenbaum, MD, FACR**
RadNet

“The reduced noise and improved image quality in the epicardial fat allows for improved depiction of the outer wall of the coronary arteries, making it easier to detect positive remodeling of high-risk plaque than with prior reconstruction methods.”

—**James Earls, MD, FSCCT**
George Washington University Hospital



Deep Learning meets an even deeper commitment.


Your dedication makes TrueFidelity truly valuable.

GE is committed to helping you succeed. That's why GE's continuing education sets the standard. Our application trainers are there to help you make the most of TrueFidelity CT Images.

It's time to get a closer look at a better way of seeing.



GE Healthcare



© 2019 General Electric Company — All rights reserved. GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. TrueFidelity CT Images is a trademark of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

Not available for sale in all regions.
JB66816XX